



Recombinant Putative G-protein coupled receptor B0244.5 (B0244.5)

Product Code	CSB-EP604645CXY-B
Abbreviation	B0244.5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10909
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MQNIFENCSYHSKYEPYFLNCTNTTNQCVLIQDVGIIQAIDFWANL CIPFTL FVIA FILN GYLSIL IPEFRKMNDTTKKQYIFVVS RGISSLSASSIMMVL RLLKMLSTSFTVYF LFFL IDDL SFY SLLGSYVGSTLLLYLATVRPIFY SIQISVRIVYKFALVNVLLAVLAVTT AIF QAAEVS DGGFFHCDVQHCQPIINIAMFVIIAT SFLIPIITLTFVLVTL CFQKSRTQSIG NF TVDNSVYKSARTR LAWTLFTFTLISLTEMIPSSFLVNL RVEDTITICVNFYQADHL FIPA IMNSFQTLAWGIALIVDPLCALLFDPRIRKVVVEHVSRLSIIIGRSFEACCHSNLN KEIQ DK
Source	E.coli
Target Names	B0244.5
Protein Names	Recommended name: Putative G-protein coupled receptor B0244.5
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.