



Recombinant Putative transporter B0252.3 (B0252.3)

Product Code	CSB-YP604646CXY
Abbreviation	B0252.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10917
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	<p>MQVRWAPSDGSLGDYTYQQDMSSSDKLSADDVLNTLDKSNRHILTCILVCGL AWSPLAFT GLCPSFVVKASENSSFIGVADEFDLTGDAWLAESTTTTFYMVGNMIGGMFIPP LADHYGR LPVFBVATVLLMAVGGMISAFSTSIMMFCIMRMIHGIFYTAAGLAGWVLGYENTP LRLRFF TSVYFGVMWVVGACFLGLLAYILPDWRYLMFCISVPNIFVALLIYMTVPESLHFL VSSQQ NEKIEAWLEKIRGPKGDISASDIVEDRDENGSSFKTLCREMWKHKMFIVYVLMV TYIWIV DTFIYFGLAFYSTNLAGNLYLNFVLMVSLVEAPAYIFSPIFMNKYGRKVLISGTHIIA GLS FLGIVLSSEAWHIHFVLLGKFAISCSFMSIYMFASEIFPTDGRNKCIGFCETLSR FGGML SPYLSHLTAVHALAPAITLSLIAVSGLLTLILPETLNTKLPSTIAETASRRQLIDD KSD SSSN</p>
Source	Yeast
Target Names	B0252.3
Protein Names	Recommended name: Putative transporter B0252.3
Expression Region	1-484
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.