



Recombinant *Oryza sativa* subsp. japonica B3 domain-containing protein Os03g0622100 (Os03g0622100, LOC_Os03g42420)

Product Code	CSB-MP604665OFG
Abbreviation	Os03g0622100, LOC_Os03g42420
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10GM4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGMLKIGKNCVCKEWQEHCYWSHMADDCKHFLTYMVGDFTESMIVPSRFA NNFNHISE VVNLKSPSGKTWSIGVAYSDTGELVLRSGWKEFVDANGVQENDCLLFRYSGV SSFDVLIF DPSGCEKASPHFVENRGFGREEKSAGAEGGGRDGDKNHGHQHQLMTPHK NSSRCRSIPS ACKRGLFSDEIEQDHREEKKEGDDEDEDEDEDEDVDKDGEDRYFPCRHGGR VTEYNLSKG DKEEISRVPVPVEPGNPVLVKVIHASHLLSSRYSTVGVVSPEFAGRYLGPAMAR EVVMERG GGGGGGDQWHVRFVRRESSRGFHGTGWRRFARDNGLLAHDVCLFELRLVD GAGAGDRLRR RPRPTMAVHVLRRVRGRFVLIR
Source	Mammalian cell
Target Names	Os03g0622100
Protein Names	Recommended name: B3 domain-containing protein Os03g0622100
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.