



Recombinant *Oryza sativa* subsp. japonica Tubby-like F-box protein 6 (TULP6)

Product Code	CSB-BP604668OFG
Abbreviation	TULP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10LG8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSFRSLIQEMRDEFSGSISRHSLRSRSHRGGGGAPRVAAVGPAEAAAMQQSC WAQLPPELL REVLVRIEESVWVPSRRDVVACAGVCRSWRGITKEIVRVPEASGKLTFFPISL KQPGPRD GTLKCFIRRNRTTQTYLYIGLTEALADDGKFLLAARKCRKPTCTDY LISLDMSD MSKGS NTYIGKLRSNFLGKFTVYDAHPPYDGAVVSKRSARVVGLNQVSPRVPAGN YPVSHISY ELNVLGARGPRRMNCIMDSIPTSAVQEGGKAPTQTEFPLSGLDSFPSISFFRS KSARIDS ATSQLSTQKEEKLVLKNKSPRWHEQLQCWCLNFRGRVTVASVKNFQLVASDE NGPTNQQE DKVILQFGKIGKDLFTMDYRYPISAFQSFACLSSFDTKIACE
Source	Baculovirus
Target Names	TULP6
Protein Names	Recommended name: Tubby-like F-box protein 6 Short name= OsTLP6 Alternative name(s): Tubby-like F-box protein 14 Short name= OsTLP14
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.