



Recombinant *Oryza sativa* subsp. japonica Probable calcium-binding protein CML27 (CML27)

Product Code	CSB-MP604670OFG
Abbreviation	CML27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10LX4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MDAAGVAAKP SLSRKPSFS RLRNGSLNAL RLRRVFDLFD RNGDGEITLD EMASALDALG LGADRAGLEA TVGGYIPAGA AGLRFGDFEA LHRALGDALF GPVEEEEPGK QGEDDDEGDM KEAFRVFDED GDGFISAAEL QAVLKKLGLP EARNLATVQE MICNVDRDCD GRVDFGEFKC MMQGITVWGA
Source	Mammalian cell
Target Names	CML27
Protein Names	Recommended name: Probable calcium-binding protein CML27 Alternative name(s): Calmodulin-like protein 27
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.