



Recombinant *Oryza sativa* subsp. japonica Ninja-family protein Os03g0214200 (Os03g0214200, LOC_Os03g11550)

Product Code	CSB-BP604678OFG
Abbreviation	Os03g0214200, LOC_Os03g11550
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10Q07
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEGFSRDLLC GIGKGGDGPR GEVRPRVDME AEEVELNLGL SLGGRFGLDR RGEKLARSSS VAAILAAPTE PSAPPSGLFR TSSLPTVAAA EAAKKQGVDE LNCRRPSGGA EAEPAAARLP ASGSPSSGSS DGEGRRLVN MDTLMRTSS LPAGIEDEWR KRKEAQLKR LEVKRKRIER RNSLTSNISK EAVGQILEEM NAGAEKVESD DDVATGNKKT GGNVNHSSDR NRCTGLPPVH RATYTQQRGS LSGIPTKHIP AMKGSADAE HNVPSAATEH RINGAIIATPP FSALAVRAVA LASRGEQLRA TGRVAARAKS MGDVERIMMQ EMPCVCTKGL PNGKRVEGFL YKYRKGEEVR IVCVCHGSFL TPAEFVKHAG GGDVANPLRH IVVNIPIPSL Y
Source	Baculovirus
Target Names	MODD
Protein Names	Recommended name: Ninja-family protein Os03g0214200
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.