



Recombinant Probable G-protein coupled receptor B0563.6 (B0563.6)

Product Code	CSB-MP604717CXY
Abbreviation	B0563.6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11082
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSLSGQLWKMVNSTSPNAAATVGNIAYCYVLPCICAIGIVGNITNLMVLASRRL RAVSYM YLRALAVADLLCMLFVLVVFVSTEYLAKNGSSINQYKLYQIQCHMLTLINWALG AGVYV VVALSLERYISIVFPMHFRTWNSPQRATRAIVIAFLIPAIFYVPYAITRYKGGKQRF DLLQ NVTIYSMDDHPIYTTFFYWQIYKWTREAILRFLPIIILTVLNIQIMIAFRKRQKMFQQ LTN KRKEQGTQKDDTLMYMLGGTVLMSLVCNIPAAINLLLIDETLKKRLDYQIFRAVA NILEI TNHASQFYVFCACSTDYRRTTFLQKFPCKTDYANRDRLRSFVRRTQSVIQKQG SVEHTTN SKVWIMFKNKFQRDSLHHSRKFSRHMPIEQDTVDIQLASGEQSTSGEMCEA DTLIKYGG TAQLCNDENNTTFL
Source	Mammalian cell
Target Names	B0563.6
Protein Names	Recommended name: Probable G-protein coupled receptor B0563.6
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.