



Recombinant Thioredoxin domain-containing protein 9 homolog (tag-170)

Product Code	CSB-BP604740CXY
Abbreviation	tag-170
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11183
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAANIQQQFG EQLLRAAQVV EEQIDQEMNK LENLEEDDLE VIRRQRMEQM KKAQKDRIEM LSHGHGKYEE VADEKEFFEA TKKSDKVVCL FYLPGNFRCK IVDKHFEILA RKHVGTRFIH VNAEKVHFLT TRLNIRVIPS IAIVVKQQTV DYIRGFDELG GKDEFTTETM ENRLARSEVL TVEKKHTAPA KKKIIRSGVE EYDNEEDW
Source	Baculovirus
Target Names	txdc-9
Protein Names	Recommended name: Thioredoxin domain-containing protein 9 homolog
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.