



Recombinant *Saccharomyces cerevisiae* Serine/threonine-protein kinase ENV7 (ENV7)

Product Code	CSB-BP604890SVG
Abbreviation	ENV7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12003
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MISIVLELFQ NLCCCRGFSD ATIRVNDKRY RIQRLLGEGG MSFVYLVQLS KNSLIIDNGI ATPELYALKK IICPSVESIS NGMREIENYK RFQSPYVIKS IDSQVMQEKD GSKTIYIVLP YYSLGSLQDS INRRLLEGTF VSEAECVRIM LGVTRGLLCL HDPASRQDNA TSRVNVDVAVS MTYSDETAML LEDTPLEMDM LSSNSAGSIA YAHRDITPSN ILFSSDGLPV IGDLGSCSQA DITIENRHQL SELQEWVNDN CTLPYTPPEL LNLKLNQVLS SKVDIWSLGC TFYTLMFGIS PFEREEQIHG ASLTYAINTG KYSPRNSRF SEGLLSVIKK CIQVDPIQRP TTSQLLNLLQ DLDT
Source	Baculovirus
Target Names	ENV7
Protein Names	Recommended name: Serine/threonine-protein kinase ENV7 EC= 2.7.11.1 Alternative name(s): Late endosome and vacuole interface protein 7
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.