



Recombinant *Saccharomyces cerevisiae* Transcriptional coactivator HFI1/ADA1 (HFI1)

Product Code	CSB-MP604900SVG
Abbreviation	HFI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12060
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSAIQSPAPK PLQPTYPAAS PASTNAYMKP GLIGSPAVSN HTEPNNGNNE TAEPQGPQR IDLGAMIEEL TSLGKESWT KYAQIISLFI LGKLSRKELS NELELVFSPS AASLEKSNTN HHHSLVRLHN QLLLGIFANS LRENPLGRNG NESSWGFGNG SNNPNNKLKR INKHNSQIEV YKKIVMSLPL NDRNRLKMIT KEAGKRGFIF CSVFQARLNN IPKIPIVTNP ESLKRVKSNN LKTPLEWSQD IMNGFNVPLA SESHSLPDTD SFYLRMVGIA REHGLVGTVD ARCVELISLA LDQYLKNIIE FTIDTVRYRR KKYSDDYDLN ESGLYKSVSE MAADKRDAKI KQLDDDKNED ECADEAKSIN NGNNSSKDDI GDISMSSITK AGEAVNEELH ENRTISLTNE DIYDLSLIFP NLVEPSGSYY ALTNLGLVND DELVDMKSNI DDLPDFLNEK PTFTPLDERN VGTRHELNLWL IKGILTED
Source	Mammalian cell
Target Names	HFI1
Protein Names	Recommended name: Transcriptional coactivator HFI1/ADA1
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.