



Recombinant *Oryza sativa* subsp. japonica U1 small nuclear ribonucleoprotein A (Os05g0154800, LOC_Os05g06280)

Product Code	CSB-MP604934OFG
Abbreviation	Os05g0154800, LOC_Os05g06280
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0DKM4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSGEVAAAVG GGAPPEENGAP PNVTIYINN L NEKIKLEELK KSLRAVFSQF GKILDVLA FK TLKHKGQAWV VFEDVASATE ALKSMQDFPF HNKPMRIQYA KTKSDIIAKA DGTFVPRERR KRNDEKPEKK QKREQHHDVS QVGLGVNAYP GVYGAPPLSQ LPFAGA QKVM MPEIIVPNNI L FVQNLPHET TPMMLQMLFC QYPGFKEVRM VEA KPGIAFV EYGDEGQATA AMNHLQGFKI TKDNQMLISY AKK
Source	Mammalian cell
Target Names	Os05g0154800
Protein Names	Recommended name: U1 small nuclear ribonucleoprotein A Short name= U1 snRNP protein A
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.