



Recombinant Schizosaccharomyces pombe DNA polymerase zeta processivity subunit (rev7)

Product Code	CSB-BP604943SXV
Abbreviation	rev7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0E7J8
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MENDTGWSVK KCIDIFGEFL LVSIHCILYA RRLYPQDLFI KARKYNTIVW QSRHPILCEY IEEVVQSCIE ELQTGSVHQV ALSIINKEQR EEERYVFSTD SIPIPDFLL EKQISTNEPF TDAYVEYMRA SLIQLLNITN GLPLIEQECT WTLRVTLKDG FPRPKQWEEW FLPPQARETD ATRQFKGITI PVRNVDIGPM MTEIWVEKYT NSD
Source	Baculovirus
Target Names	rev7
Protein Names	Recommended name: DNA polymerase zeta processivity subunit Alternative name(s): Revertibility protein 7
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.