



Recombinant *Oryza sativa* subsp. *japonica*

Probable protein phosphatase 2C 67

(Os09g0314400, LOC_Os09g14540)

Product Code	CSB-EP605203OFG-B
Abbreviation	Os09g0314400, LOC_Os09g14540
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0J2R1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MAHQKREATS DNGGGDEEWA SKRPKVVGAA AEKEHILTSD ASHETNGDEA QGGDASRKEN TVSTNPCVSD EKAATNSNVS SGHGVILTSV EADAAEDKGC RHTMEDAWVL LPDASMESPG NLRCAHFAYI DGHGGRLAAE YAQKHLHQNV IAAGLPRELM DVKAAKKAI EGFRRTDECL LQUESTKGNWQ DGATAVCVWV LGQTVVVANA GDAKAVLARS TSADGEGAVD DAKSQLKAI VLTREHKAIFP QERARIQKAG GSVGPNGLRLQ GRIEVSRLG DRQFKKVGLI ATPDVHSFEV TRKDHFIILG CDGLWGVFGP GDAVEFVQNN LKETSSATLA VRRLVKEAVR ERRCKDNCTA VLIVFKH
Source	E.coli
Target Names	Os09g0314400
Protein Names	Recommended name: Probable protein phosphatase 2C 67 Short name= OsPP2C67 EC= 3.1.3.16
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.