



Recombinant *Oryza sativa* subsp. *japonica* Zinc finger CCCH domain-containing protein 28 (Os04g0438700, LOC_Os04g35800)

Product Code	CSB-EP605210OFG
Abbreviation	Os04g0438700, LOC_Os04g35800
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JD07
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MASAETPNPD AEIPNTDAAA AADPAAAAPA AAATDPAAAG SPSPLPPRK RRLSPTPSPT RRSSRSRSRS PRRGRSRSRS RSRSRGRSAS PRYPDGKRRR HNDLNVEVCR DFLRDRCARA DIECKYAHPH PTVAVDRDSK VTACADSLRN NCFRGRTRY YHPPPHIQES LLRSIGVEDP KVKMQVCRDF TRGRCSRSAN ECRFLHHSPL EDCAIVCQDF LRGRCDRKSC RYSHVMAHPM PPPMRDIPMQ YPDMVYMPPP APLGVPMMP PPSAPAAFSG NNYGVEVCRD YLKNMCNRES CRFAHPDLNN EVMNTQVEVC RDFKRGEENR PACRFYHPPA SSNSIG
Source	E.coli
Target Names	Os04g0438700
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 28 Short name= OsC3H28
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.