



# Recombinant *Cupriavidus necator* Fructose-1,6-bisphosphatase class 1 1 (fbp1)

<b>Product Code</b>	CSB-BP605255DZQ
<b>Abbreviation</b>	fbp1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0KCY0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Cupriavidus necator</i> (strain ATCC 17699 / H16 / DSM 428 / Stanier 337) ( <i>Ralstonia eutropha</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTRISLTRYL VEEQRKHNTI QPELRLLIEV VARACKAISN AVSKGALAGV LGSAGTGNVQ GETQQKLDVI ANEVLLDANE WGGHLAAMAS EEMESFYEIP NRYPKGEYLL MFDPLDGSSN IDVNVSIGTI FSVLHMPKPG QTVTEADFMQ PGTHQVAAGY AVYGPQTTLV LTVGNGVHMF TLDREAGSFV LTHSNVTIPD DTKEFAINMS NMRHWAPPVR RYIDECLAGE EGPRGKNFNM RWWASMVADV HRILTRGGVF MYPWDKREPE KPGKLRLMYE ANPMAMLVEQ AGGAATNGHQ RIMDVQPEKL HQRVSVILGS KNEVERVTRY HLEAEDKA
<b>Source</b>	Baculovirus
<b>Target Names</b>	fbp1
<b>Protein Names</b>	Recommended name: Fructose-1,6-bisphosphatase class 1 1 Short name=FBPase class 1 1 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6-bisphosphate 1-phosphohydrolase class 1 1
<b>Expression Region</b>	1-338
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.