



Recombinant Gorilla gorilla gorilla NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 5 (NDUFA5)

Product Code	CSB-EP605270GGZ-B
Abbreviation	NDUFA5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0MQA1
Product Type	Recombinant Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Purity	>85% (SDS-PAGE)
Sequence	AGVLKKTG LVGLAVCNTP HERLRILYTK ILDVLEEIPK NAAYRKYTEQ ITNEKLAMVK AEPDVKKLED QLQGGQLEEV ILQAEHELNL ARKMREWKLW EPLVEEPPAD QWKWPI
Source	E.coli
Target Names	NDUFA5
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 5 Alternative name(s): Complex I subunit B13 Complex I-13kD-B Short name= CI-13kD-B NADH-ubiquinone oxidoreductase 13 kDa-B subunit
Expression Region	2-116
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.