



Recombinant Escherichia coli O6:K15:H31 50S ribosomal protein L11 (rpLK)

Product Code	CSB-EP605565EGY
Abbreviation	rpLK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TA82
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	AKKVQAYVK LQVAAGMANP SPPVGPALGQ QGVNIMEFCK AFNAKTDSIE KGLPIPVVIT VYADRSFTFV TKTPPAAVLL KKAAGIKSGS GKPNKDKVVGK ISRAQLQEIA QTKAADMTGA DIEAMTRSIE GTARSMGLVV ED
Source	E.coli
Target Names	rpLK
Protein Names	Recommended name: 50S ribosomal protein L11
Expression Region	2-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.