



# Recombinant Escherichia coli O6:K15:H31 sn-glycerol-3-phosphate import ATP-binding protein UgpC (ugpC)

<b>Product Code</b>	CSB-EP605581EGY
<b>Abbreviation</b>	ugpC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TC10
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGLKLQAVT KSWDGKTQVI KPLTLDVADG EFIVMVGPSG CGKSTLLRMV AGLERVTTGD IWIDRKRVT MEPKDRGIAM VFQNYALYPH MSVEENMAWG LKIRGMGKQQ IAERVKEAAR ILELDGLLKR RPRELSGGQR QRVAMGRAIV REPAVFLFDE PLSNLDAKLR VQMRLELQQL HRRLKTTSLY VTHDQVEAMT LAQRVMVMNG GVAEQIGTPV EYEEKPASLF VASFIGSPAM NLLAGRVNNE GTHFELDGGI TLLLNGGYRQ YAGRKMTLGI RPEHIALSSR AEGGVPLVMD TLEILGADNL AHGRWGEQKL VVRLAHQERP TAGSTLWLHL PENQLHLFDG ETGQRV
<b>Source</b>	E.coli
<b>Target Names</b>	ugpC
<b>Protein Names</b>	Recommended name: sn-glycerol-3-phosphate import ATP-binding protein UgpC EC= 3.6.3.20
<b>Expression Region</b>	1-356
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.