



Recombinant Escherichia coli O6:K15:H31 Putrescine aminotransferase (patA)

Product Code	CSB-MP605595EGY
Abbreviation	patA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TD34
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNRLPSSASA LAYSAHALNL IEKRTL DHEE MKALNREVIE YFKEHVNPGF LEYRKSVTAG GDYGAVEWQA GGLNTLVDTQ GQEFIDCLGG FGIFNVGHRN PVVSAVQNNQ LAKQPLHSQE LLDPLRAMLA KTLAALTPGK LKYSFFCNSG TESVEAALKL AKAYQSPRGK FTFIATSGAF HGKSLGALSA TAKSTFRKPF MPLLPGFRHV PFGNIEAMRT ALSECKKTGD DVA AVILEPI QGEGGVILPP PGYLTAVRKL CDEFGALMIL DEVQTMGRT GKMFACHEHEN VQPDI LCLAK ALGGGVMPIG ATIAATEEVFS VLFDNPFLHT TTFGGNPLAC AAALATINVL LEQNLPAQAE QKGDMLLDGF RQLAREY PDL VQEARGKGML MAIEFVDNEI GYNFASEMFR QRVLVAGTLN NAKTIRIEPP LTLTIEQCEQ VIKAARKALA AMRVSVEEA
Source	Mammalian cell
Target Names	patA
Protein Names	Recommended name: Putrescine aminotransferase EC= 2.6.1.82 Alternative name(s): Putrescine--2-oxoglutaric acid transaminase Short name= PAT Short name= PATase
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.