



Recombinant Escherichia coli O6:K15:H31 Betaine aldehyde dehydrogenase (betB)

Product Code	CSB-MP605650EGY
Abbreviation	betB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TKW0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSRMAEQQLY IHGGYTSATS GHTFETINPA NGNVLATVQA AGREDVDRAV KSAQQGQKIW AAMTAMERSR ILRRAVDILR ERNDELAKLE TLDTGKAYSE TSTVDIVTGA DVLEYYAGLI PSLEGSQIPL RETSFVYTRR EPLGVVAGIG AWNYPHQIAL WKSAPALAAG NAMIFKPSEV TPLTALKLAE IYSEAGLPDG VFNVLPGVGA ETGQYLTEHP GIAKVSFTGG VASGKKVMAN SAASSLKEVT MELGGKSPLI VFDDADLELA ADIAMMANFF SSGQVCTNGT RVFVPAKCKA AFEQKILARV ERIRADDVFD PQTNFGPLVS FPHRDNVLRY IAKGKEEGAR VLCGGDVLKG DGFNGAWVA PTVFTDCSDD MTIVREEIFG PVMSILYIES EDEVIRRAND TDYGLAAGIV TADLNRAHRV IHQLEAGICW INTWGESPAAE MPVGGYKHSG IGRENGVMTL QSYTQVKSIQ VEMAKFQSIF
Source	Mammalian cell
Target Names	betB
Protein Names	Recommended name: Betaine aldehyde dehydrogenase Short name= BADH EC= 1.2.1.8
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.