



Recombinant Bovine Transcription factor jun-B (JUNB)

Product Code	CSB-BP605746BO
Abbreviation	JUNB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0VBZ5
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MCTKMEQPFY HDDSYAAAGY GRTPGGLSLH DYKLLKPSLA LNLSDPYRNL KAPGARGPGP EGNNGGSYFS SQGSDTGASL KLASSELERL IVPNSNGVIT TTPTPPGQYF YPRGGGSGGG AGGAGGGVTE EQEGFADGFV KALDDLHKMN HVTTPNVSLG ASGGPPAGPG GUYAGPEPPP VYTNLSSYSP ASAPSGGAGA AVGTGSSYPT ATISYLPHAP PFAGGHPAQL GLGRGASAFK EEPQTVPEAR SRDATPPVSP INMEDQERIK VERKRLRNRL AATKCRKRKL ERARLEDKV KTLKAENAGL SSTAGLLREQ VAQLKQKVM T HVSNGCQLLL GVKGFHAF
Source	Baculovirus
Target Names	JUNB
Protein Names	Recommended name: Transcription factor jun-B
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.