



Recombinant Human Protein INCA1 (INCA1)

Product Code	CSB-MP605764HU
Abbreviation	INCA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0VD86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQVQDDGVNL IPFAKCSRVV SRSPPPRLPS QSLRPMPQRY GDVFWKNLNQ RPTPTWLEEQ HIPPMRLRATG CSQLGLYPPE QLPPPEMLWR RKKRRPCLEG MQQQGLGGVP ARVRAVITYHL EDLRRRQSII NELKKAQWGS SGAASEPVVL GEEGCGFPST NEYPDLEER ATYPQEEDRF LTPGRAQLLW SPWSPLDQEE ACASRQLHSL ASFSTVTARR NPLHNPWGME LAASEE
Source	Mammalian cell
Target Names	INCA1
Protein Names	Recommended name: Protein INCA1 Alternative name(s): Inhibitor of CDK interacting with cyclin A1
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.