



Recombinant Arabidopsis thaliana Probable S-acyltransferase At3g26935 (At3g26935)

Product Code	CSB-EP605883DOA
Abbreviation	At3g26935
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0WQK2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MYVPPPPQRSDSGNSGDLRVYQTWKGSNIFFLQGRFVFGPDVRSALTICLIA VPVTIFC IFVARKLIDDFSDSWGVSIVAVAVVFTIYDLILLLLTSGRDPGIIPRNAHPPEPETL DGN MDAGAGQTPQLRLPRIKEVQLNGITFKVKYCDTCMLYRPPRCSHCSICNNCVE RFDHHCP WVGQCIGMRNYRFFFMFVFSTLLLCIYVFAFCWVYIRKIMSEHTTTWKAMLK TPASIVL IITYFISMWFVGGTLVFHLYLISTNQTTYENFRYRYDRRSNPHNKGVVNNFKET FFSTIP PSKNDFRAMVQREPPLPPRSVAGGFMSPNMGKANDEIEMGRKAVWADMGP AMSDHGDGKH GNNERLHVKDGELGELSPDIRATVDEQSDRPSMHPRRSSWGRKSGSWDMS PEVMALAARV GEQNQNGGGSSSGSGLVTENRPT
Source	E.coli
Target Names	PAT07
Protein Names	Recommended name: Probable S-acyltransferase At3g26935 EC= 2.3.1.- Alternative name(s): Probable palmitoyltransferase At3g26935 Zinc finger DHHC domain-containing protein At3g26935
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.