



Recombinant Arabidopsis thaliana Probable ADP-ribosylation factor GTPase-activating protein AGD15 (AGD15)

Product Code	CSB-BP605884DOA
Abbreviation	AGD15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0WQQ1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNRKASVSKE LNAKHSKILE ALLKHPDNRE CADCRSKAPR WASVNLGIFI CMQCSGIHRS LGVHISQVRS ITLDTWLDPDQ VAFMKSTGNA KGNEYWESEL PQHFERSSSD TFIRAKYSEK RWWSPGAIQP APIVSQLSCK VSHLVESGYK PETPKKARTL SLDEEILLHH VLQVTPPETR TRAGSVDMEKE NVYVVPLPEF KKPNQKNENF SSEVNQNRRT TIAPPSSWAT FD
Source	Baculovirus
Target Names	AGD15
Protein Names	Recommended name: Probable ADP-ribosylation factor GTPase-activating protein AGD15 Short name= ARF GAP AGD15 Alternative name(s): Protein ARF-GAP DOMAIN 15 Short name= AtAGD15
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.