



Recombinant Arabidopsis thaliana F-box protein SKIP17 (SKIP17)

Product Code	CSB-MP605885DOA
Abbreviation	SKIP17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0WRC9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDPPTQSHNS VRTVPGSNHS HTDLISSFL SFPDSSPISI SNSFDRVLDR ALASASADES VQDRLVDRTL ELASLLLDST RRCFRKRASV HNSNSWSLPP ELTIKVFMSL DTKSMMQAAV CCTMFNKCAM DRLCYSHIDL TTSARYADKG VVSTMINRAG KELRSLKLGR VVRTAGSDSA APLLSGSCLS PLAYNHGFLG SRLRSLRLYN LRPIKYRSLC DALSVCPNIT DLRIVGLYNL TEELFNSLTK KCRLIEHLFL ETYGYPRLE SKAGSSLVEF VTNCPNLTSL TLIRFGLTDD WARNLAESCR KLKYLNLSRS PTIKGRFLRE LGLSCKENLL KTLILRSCPK LQEKEVLEFC NSLLTGNFKS IRQIDVSSNS GLASSDRGKR CNKPNFPLER LKEERSDVTF VADFPSTSSG KRYGVCDEEE LRLIEMMEAE DDEVDEEDDS DDDTDDVSDE DESENDDDMG MGFVDVYLL
Source	Mammalian cell
Target Names	SKIP17
Protein Names	Recommended name: F-box protein SKIP17 Alternative name(s): SKP1-interacting partner 17
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.