



Recombinant Schizosaccharomyces pombe Nucleoporin seh1 (seh1)

Product Code	CSB-MP605935SXV
Abbreviation	seh1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10099
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDDISATTIQ TNHQDLVNDV TYDFYGRMV SCSADQRVKV YDFNDDTETW AITSEWRAGD ASLMRVAAH PSFGQVLAVC SLDRGVRIYE EQKKNFESKT WVEVAKLMDA RSAVLDISFC PFQHGCKLAA VSADATLRIY EAMEPGNLTY WTLMNEIALM PSPPSRNEQP AFCVNWCP SR WREQYIAVGC MNDAYIYKQN SHGKWKVVAE LPGHTDLIRD ICWAPSMGSS YYLIATAACKD GNVRIKQVET LCEEVFQEEE DAGNSMTEDS NFNLNSLKVE LIGEYDNHKC QVWRCRFNVT GTILSSSGDD GCVRLWKASY ANLFKCSVV SLEKKPEKL
Source	Mammalian cell
Target Names	seh1
Protein Names	Recommended name: Nucleoporin seh1 Alternative name(s): Nuclear pore protein seh1
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.