



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C3H8.04 (SPAC3H8.04)

Product Code	CSB-MP605943SXV
Abbreviation	SPAC3H8.04
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10140
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MPACYCATST SLENNELLYK NIRNLFSTSR SFPIEQEWMN LKSISQMKDF FSNFPGNSKA NNHFLCNSPL KFEIFNNEKS VKPSNGPHLF TRCSCRCNKL LSGDYITQQQ PIDASAEHSL NAKGINTYSL KNCFQRSGFI GSYEESLFSG HMPYCSSLVVF EFSIEIGVIS FCRCKPSLVF PPHLKINFVA YSLVGNVDNV QFPYIGRFRL RSQKSDKVMN KGYPFGYRIP SVGQLQLILR QTNGLVIKVF LVPYNVSSMV DCSKTWIRQK HYLQQLDDKS GKILSHLKFQ LQLQIICTSA GHHYLYDSQR IIFVQQSLGG LYGNTKIVNE TLLSESCR
Source	Mammalian cell
Target Names	SPAC3H8.04
Protein Names	Recommended name: Uncharacterized protein C3H8.04
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.