



Recombinant Schizosaccharomyces pombe Deubiquitination-protection protein dph1 (dph1)

Product Code	CSB-MP605947SXV
Abbreviation	dph1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10169
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTNISLTIKA ANDQKYAVTV DSESSVLALK EAIAPVADIE KERQRLIYAG RVLKDEESLK TYKIQDGHSI HLVKTLGQNP AAAATNVSDR TQQVPTNIQA GGGANNPLAN LTSARYAGFN IPMPASMFG PNPENPVPPS TEELANMLSN PMVQSSINEM FSNPQMLDMI INSSPHLRNA PPYVRQMMQS PEFRRAMTDP DTMRQMAQLH QQMGAAGIDP MSLMGGGLGG AGLGGLGGAG LGGFGGANNA TAGIAGAAPV DQTAAANTIQLLLNNLGGAG FGAGLGDAGL GAGLGGAASP PAPAQDTRPP EERYAEQLSQ LNEMGFVDFE RNVQALRRSG GNVQGAIESL LSDL
Source	Mammalian cell
Target Names	dph1
Protein Names	Recommended name: Deubiquitination-protection protein dph1
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.