



Recombinant Schizosaccharomyces pombe RNA polymerase I-specific transcription initiation factor RRN9 homolog (acr1)

Product Code	CSB-BP605954SXV
Abbreviation	acr1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10204
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNDSEEELEQ QGSQEDISSS QDTSDEKEMR KQDALYAQLI EYREIVSELR RSEMKNLCIQ LYCSYKAGIT RRNELQSDLD GYERRIPRL RAVWPTIHSI PDESFFDGTS ASRFRQEAY LAVRLASKKL HKLGRMASSE ELPNPDGLRY LVERLQIKLE RLFDSVDTFR QEQGSKSMQS RVKPMDWKIL YGLARVTNGS FSSSQKADVL ALGNGTRACQ SLFKEMSLQP SIDIEKVSHS NPEISIEFED
Source	Baculovirus
Target Names	acr1
Protein Names	Recommended name: RNA polymerase I-specific transcription initiation factor RRN9 homolog Alternative name(s): Accumulation of condensin at rDNA protein 1
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.