



Recombinant *Schizosaccharomyces pombe* Uncharacterized J domain-containing protein C4H3.01 (SPAC4H3.01)

Product Code	CSB-BP605955SXV
Abbreviation	SPAC4H3.01
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10209
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MATPVDTEYY DLLGISTDAT AVDIKKAYRK LAVKYHPDKN PDDPQGASEK FQKISEAYQV LGDEKLR SQY DQFGKEKAVP EQGFTDAYDF FTNLFGGAPF REWVGELSFV KEMFREEDSA VEQQQMNDKQ QLLLESSEPT PTIKQQFNDR KKNAQIRERE ALAKREQEMI EDRRQRIKEV TENLEKRLDD WIAKATTEEG LNALREKYTQ EANTLRIESF GVEILHAIGE VYTQKGRTVL KSSKFGIGGF WSRMKEKGKI ARATWDTVSA AMDAKLSIDQ MQKLEDKGED QASAEERAKL ELDITGKILR ASWCGARYDI QGVLREACSN LLKKRVPTTEL RLKRAHALLE IGTIFSNVEA DPDDPNRIFE NLILENKKKR KKGSEAKPPK AT
Source	Baculovirus
Target Names	SPAC4H3.01
Protein Names	Recommended name: Uncharacterized J domain-containing protein C4H3.01
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.