



Recombinant Schizosaccharomyces pombe 60S ribosomal protein L9-A (rpl9a)

Product Code	CSB-MP605959SXV
Abbreviation	rpl9a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10232
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGRDIYKDET LTIPKGVTV D IKARNVTVTG PRGTLKQNL R HVDIEMKKQG NTIKFIVWHG SRKHNACIRS VYSIINNMII GVTQGFYK M RLVYAHFPIN INLTENGT VV EIRNFLGERI TRVIKCLPGV TVSISSAVKD EIILEGNSLE NVSQSAANIK QICNVRNKDI RKFLDGIYVS ERGNIEELE E
Source	Mammalian cell
Target Names	rpl901
Protein Names	Recommended name: 60S ribosomal protein L9-A
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.