



Recombinant Human D site-binding protein (DBP)

Product Code	CSB-BP606009HU
Abbreviation	DBP
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q10586
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MARPVSDRTP APLLLGGPAG TPPGGGALLG LRSLLQGTSK PKEPASCLLK EKERKAALPA ATTPGPGLET AGPADAPAGA VVGGGSPRGR PGPVPAPGLL APLLWERTLP FGDVEYVDLD AFLLEHGLPP SPPPPGGPSP EPSPARTPAP SPGPGSCGSA SPRSSPGHAP ARAALGTASG HRAGLTSRDT PSPVDPDTVE VLMTFEPDPA DLALSSIPGH ETFDPRRHRF SEEELKPQPI MKKARKIQVP EEQKDEKYWS RRYKNNEAAK RSRDARLKE NQISVRAAFL EKENALLRQE VVAVRQELSH YRAVLSRYQA QHGAL
Source	Baculovirus
Target Names	DBP
Protein Names	Recommended name: D site-binding protein Alternative name(s): Albumin D box-binding protein Albumin D-element-binding protein Tax-responsive enhancer element-binding protein 302 Short name= TaxREB302
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>