



Recombinant G2/mitotic-specific cyclin-B3 (cyb-3)

Product Code	CSB-EP606013CXY-B
Abbreviation	cyb-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10654
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MMLRSQAKNV DLTSQADSRH QQKRKQAEQL DALKNPSEPA AKKQHSKGLT ELRAHISGFK IDSAKRDPLG KSRTSRRDVE NLPPQKSRYV DPCPHYDYDL EEAGNPDSIS DYAQGIFDYY RHREVFHFRVR KYLHKHPEVD VKTRAILIDW MVEIQETFEL NHETLYNAVK LTDMYLCKTK NVDKNTIQKL ACVAIFIAAK YDERSPPPLVD DLIYLSGDRF SRDELLAMER ELFATVGYDL GSPLSYRYLR RFGRVCRVDM KTLTMGRFIL ETSLMVYEYA MVSQSRLAAA AFVLAMRMLD KNNEYEWNPV LEKYSGFTGE EVMPLVEHMN HILHFSKDKW AQLTSVRQKY SHEVFFHVAS IPMLPDTLKV VDSHTYAPVP MLSYP
Source	E.coli
Target Names	cyb-3
Protein Names	Recommended name: G2/mitotic-specific cyclin-B3
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.