



Recombinant Uncharacterized protein F22D3.4 (F22D3.4)

Product Code	CSB-EP606040CXY
Abbreviation	F22D3.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10919
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNLYLLVLSSTFLCLAAEEIYKTNGNTTISENSNILNVLPNSDIWKTGSITHDKKD YLEA LFSKSYKRECEHETRLEPVYFYPGENVKIKCYMCNLALVFNGNPKGWARVK NIENFLEK GEKFGSKDGPAAEMVENS DIVEDQNAQSSARSLGKPYYYEEDGFLVIENANV YSQGVYFC FDEDSVASQRFFYILIPILPVFHIDSKNNSKIMELTNNCSDSVRKFPNHYYMNS NEKVPF DDQDTCFARNGMDDWKCYDFEKIKRFEKCNPATSCVKKFVHPTLPNDIPISI NLQWSPW SECNDRTQTRTGQCFIRLERQLDDSDGLPENWSWLRKLNTMEKHLAFRDGV PLFNGLLAS WLYEVESLESCLDTELGPKSKLSGYDNFFGNVLTSLFSNKTISGIKMSDNAFKY CIRTKT LDETNDALIGTYTEDKRTC
Source	E.coli
Target Names	F22D3.4
Protein Names	Recommended name: Uncharacterized protein F22D3.4
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.