



Recombinant Putative acid phosphatase 5 (pho-5)

Product Code	CSB-EP606045CXY
Abbreviation	pho-5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10944
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	INAVVYK EVPIQANTDT LEYVHTVWRH GDRTPAELLF PDDITKWPEG LGELTEQGAA QQYRLGQWLK RRYGSQLGK FNRNAIYIRS SDYNRTLMSA QANMAGLFPP KYPIAGGLMW QPIPVHTISK PTDKELYEEA SCPTAEIEMN AQWKSTKANG IRKKFARELS FFSQKLNLPN MELKATWRIF DNLFCQKQNN ITWPSWMNSS IFERVDQLYN EVSQLEFHTD TLRRLRGGTL LEEIFHRFSD KASGSLGKEA KFYAYSAHDS TIAALLATLG VFYDIYPKYA TCLLIEMHKL ANETRLIRVF HKNETDIDRL IEYSIPGCDD PCTLQKLGDD LKKYFPEDWE AECGLKTSFQ FIYLVIIISIL VISTVCSCCTM LFVEKQKRKI LRFPVDGLRD DTAPMLGGDD SD
Source	E.coli
Target Names	pho-5
Protein Names	Recommended name: Putative acid phosphatase 5 EC= 3.1.3.2
Expression Region	14-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.