



# Recombinant Phosphatidylserine decarboxylase proenzyme (psd-1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP606046CXY   |
| <b>Abbreviation</b>      | psd-1   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q10949  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Caenorhabditis elegans  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MMPLFNVLRS ARMLPAVSKK VVSPPMMLRS VRELTNQSKN VYATKEVIIG ASQKKKRSWV KWLSVSTLII GGASYVGYLF TPDWREIVDS KHYYSNWKIR VYLSLPFNTA SRVIGGLANQ EIPVWLREHL LGGFARMYDC RMDDCVDPDF KNYPFAAFF NRKLKESTRP ISASPLVSPA DGTVLHFGKV EDNKIEYVKG HDYDVKFLG DVDLPQKDEL DLYQVVIYLA PGDYHAFHSP ARWVANQCRH VPGLLSVRP TLLSHVPHLF CLNERVVLNG SWRHGFFSMS AVAATNVGDI VVDAEPSLRT NIVRRKTQKI MNTETEIHAP YVSGERVGEF RLG                      |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | psd-1   |
| <b>Protein Names</b>     | Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain   |
| <b>Expression Region</b> | 1-343   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |