



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-SDS-like (Os03g0225200, LOC_Os03g12414)

Product Code	CSB-EP606068OFG-B
Abbreviation	Os03g0225200, LOC_Os03g12414
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10PQ9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MPPTMLASVP TRPRSHPFRR RRGAAAAAPP LLPDQIAAAA AAAAKRPAES STSASSCFHS EVISATSTTC PTSLAAAQRP EKRPRYQDVD EEQPAASECS EIIGGARPRA AEVEVSESSC LASVLESYLA CPEQLANDAE TTAYSSARED LTLSETEEEE EEEVRSRSGPC ICTDCSFSPL HESSSSSDDD NAVPSPTFSL FLALAEQFVP FTHPKTPTAT DVALQAGEGK RFEDLDNEVS YERFRRRERR GVVARDYIEV YSSMLGSYGR AVVEQRVVMV NWIMEHSQAM KLQPETVFMG IGLMDRFLTR GYVKGSRNLQ LLGIACTTLA TRIEENQPYN CILQKAFKVG INTYSRSEVV AMEWLVQEV LDFQCFVTTTH HFLWFYKAA NADDRVEDLA KYLALLSLLD HKHLSFWPST VAAAVVALAC LATNNESSCH LVMETHMRTK NDDLPECLMS LEWLTNYAS
Source	E.coli
Target Names	SDS
Protein Names	Recommended name: Cyclin-SDS-like Alternative name(s): Protein SOLO DANCERS-like
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.