



# Recombinant *Oryza sativa* subsp. *japonica* Cyclin-SDS-like (Os03g0225200, LOC\_Os03g12414)

<b>Product Code</b>	CSB-EP606068OFG
<b>Abbreviation</b>	Os03g0225200, LOC_Os03g12414
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q10PQ9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPPTMLASVP TRPRSHPFRR RRGAAAAAPP LLPDQIAAAA AAAAKRPAES STSASSCFHS EVISATSTTC PTSLAAAQRP EKRPRYQDVD EEQPAASECS EIIGGARPRA AEVEVSESSC LASVLESYLA CPEQLANDAE TTAYSSARED LTLSETEEEEE EEEVVRSGPC ICTDCSFSPL HESSSSSDDD NAVPSPTFSL FLALAEQFVP FTHPKTPTAT DVALQAGEGK RFEDLDNEVS YERFRRRERR GVVARDYIEV YSSMLGSYGR AVVEQRVVMV NWIMEHSQAM KLQPETVFMG IGLMDRFLTR GYVKGSRNLQ LLGIACTTLA TRIEENQPIN CILQKAFKVG INTYSRSEVV AMEWLVEVL DFQCFVTTTH HFLWFYKAA NADDRVEDLA KYLALLSLLD HKHLSFWPST VAAAVVALAC LATNNESSCH LVMETHMRTK NDDLPECLMS LEWLTYAS
<b>Source</b>	E.coli
<b>Target Names</b>	SDS
<b>Protein Names</b>	Recommended name: Cyclin-SDS-like Alternative name(s): Protein SOLO DANCERS-like
<b>Expression Region</b>	1-469
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.