



Recombinant Putative serine/threonine-protein kinase C05D10.2 (C05D10.2)

Product Code	CSB-EP606129CXY-B
Abbreviation	C05D10.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11179
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTDDVDTHIHEKFDLQKRLGKGAYGIVWKAYDKRSRETVALKKIFDAFRNPTD SQRTFRE VMFLQEFGKHPNVIKLYNIFRADNDRDIYLAFEFMEADLHNVIKKGSIKDVHKQ YIMCQ LFRAIRFLHSGNVLHRDLKPSNVLLDADCRVKLADFGGLARSLSSLEDYPEGQK MPDLTEY VATRWYRSPEILLAAKRYTKGVDMWSLGCILAEMLIGRALFPGSSTINQIERIMN TIAKP SRDIASIGSHYAASVLEKMPQRPRKPLDLITQSQTAAIDMVQRLLIFAPQKRL TVEQC LVHPYVVQFHNPSEEPVLNYEVYPPPLPDHIQLSIDDYRDRLYEMIDEKKASFKRI QHEKI RPYGEDKSRAPIAQAECSDTDYDTARSLQRTTSMKNNSSSHDSSSGTLRER AASAESRT SKDSNGEMRNGNGNTPSSIKQRRRSVERARLFANIKPSKILHPHKLISNY
Source	E.coli
Target Names	mapk-15
Protein Names	Recommended name: Putative serine/threonine-protein kinase C05D10.2 EC=2.7.11.24
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.