



# Recombinant *Oryza sativa* subsp. *japonica* Guanine nucleotide-binding protein alpha-1 subunit (GPA1)

<b>Product Code</b>	CSB-EP606324OFG-B
<b>Abbreviation</b>	GPA1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0DJ33
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GSSCSRSHS LSEAETTKNA KSADIDRRIL QETKAEQHIH KLLLLGAGES GKSTIFKQIK LLFQTGFDEA ELRSYTSVIH ANVYQTIKIL YEGAKELSQV ESDSSKYVIS PDNQEIGEKL SDIDGRLDYP LLNKELVLDV KRLWQDPAIQ ETYLGRSILQ LPDCAQYFME NLDRLAEAGY VPTKEDVLYA RVRTNGVVQI QFSPVGENKR GGEVYRLYDV GGQRNERRKW IHLFEGVNAV IFCAAISEYD QMLFEDETKN RMMETKELFD WVLKQRCFEK TSFILFLNKF DIFEKKIQKV PLSVCEWFKD YQPIAPGKQE VEHAYEFVKK KFEELYFQSS KPDRVDRVFK IYRTTALDQK LVKKTFLKID ESMRRSREGT
<b>Source</b>	E.coli
<b>Target Names</b>	GPA1
<b>Protein Names</b>	Recommended name: Guanine nucleotide-binding protein alpha-1 subunit Short name= GP-alpha-1 Alternative name(s): Protein Dwarf1
<b>Expression Region</b>	2-380
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.