



Recombinant Danio rerio Phosphotriesterase-related protein (pter)

Product Code	CSB-MP606681DIL
Abbreviation	pter
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0P3Z2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSELRGKVQT VLGQIEPDQL GRTMTHEHLT MSFECSHVPP APGDEGLSLA PIEMKHLHWL QQNPYSHHEN LLLNQELEAV KEELLCYRKA GGGTIVENTT TGINRNLPAL KQLAKETGVH VIAGAGYYVD VTHSEETKRM TVEKLTDVIV SEVLHGADGT DIRCGVIGEI GTSWPITese KKV LQATAHA QTRLGCPVII HPGRDNRAPV EVIRILQEAG GDISKTVM SH LDRSIYDHGE LLEFARMGSY LEYDLFGTEV LNYQFNCNVD MP SDSQRVQS LKFLIQEGYE DRILIAHDIH TKHRLTKYGG HGF SHILKNI VPKMLSRGIT QNQVDKILIE NPKRWLTFK
Source	Mammalian cell
Target Names	pter
Protein Names	Recommended name: Phosphotriesterase-related protein EC= 3.1.-.- Alternative name(s): Parathion hydrolase-related protein
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.