



Recombinant Arabidopsis thaliana Aldo-keto reductase family 4 member C9 (AKR4C9)

Product Code	CSB-MP606708DOA
Abbreviation	AKR4C9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0PGJ6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MANAITFFKL NTGAKFPSVG LGTWQASPLG VGDVAAAVK IGYRHIDCAQ IYGNEKEIGA VLKCLFEDRV VKREDLFITS KLWCTDHDPQ DVPEALNRTL KDLQLEYVDL YLIHWPARIK KGSVGIKPEN LLPVDIPSTW KAMEALYDSG KARAIGVSNF STKKLADLLE LARVPPAVNQ VECHPSWRQT KLQEFCKSKG VHLSAYSPLG SPGTTWLKSD VLKNPILNMV AEKLGKSPAQ VALRWGLQMG HSVLPKSTNE GRIKENFNVF DWSIPDYMFA KFAEIEQARL VTGSFLVHET LSPYKSIEEL WDGEI
Source	Mammalian cell
Target Names	AKR4C9
Protein Names	Recommended name: Aldo-keto reductase family 4 member C9 EC= 1.1.1.-
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.