



Recombinant Escherichia coli O6:K15:H31 Aspartate carbamoyltransferase (pyrB)

Product Code	CSB-MP606950EGY
Abbreviation	pyrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0T9E3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MANPLYQKHI ISINDLSRDD LNLVLATAAK LKANPQPELL KHKVIASCFF EASTRTRLSF ETSMHRLGAS VVGFSDSANT SLGKKGETLA DTISVISTYV DAIVMRHPQE GAARLATEFS GNVPVLNAGD GSNQHPTQTL LDLFTIQETQ GRLDNLHVAM VGDLKYGRTV HSLTQALAKF DGNRFYFIAP DALAMPQYIL DMLDEKGIW SLHSSIEEVM AEVDILYMTR VQKERLDPSE YANVKAQFVL RASDLHNAKA NMKVLHPLPR VDEIATDVDK TPHAWYFQQA GNGIFARQAL LALVLNRDLV L
Source	Mammalian cell
Target Names	pyrB
Protein Names	Recommended name: Aspartate carbamoyltransferase EC= 2.1.3.2 Alternative name(s): Aspartate transcarbamylase Short name= ATCase
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.