



Recombinant Escherichia coli O6:K15:H31 Cytochrome c-552 (nrfA)

Product Code	CSB-EP606953EGY
Abbreviation	nrfA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0T9X8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	EQTA APATPVTVEA KNETFAPQHP DQYLSWKATS EQSERVDALA EDPRLVILWA GYPFSRDYNK PRGHAFVTD VRETLRTGAP KNAEDGPLPM ACWSCKSPDV ARLIQKDGED GYFHGKWARG GPEIVNNLGC ADCHNTASPE FAKGKPELTL SRPYAARAME AIGKPFKAG RFDQQSMVCG QCHVEYYFDG KNKAVKFPWD DGMKVENMEQ YYDKIAFSDW TNSLSKTPML KAQHPEYETW TAGIHGKNNV TCIDCHMPKV QNAEGKLYTD HKIGNPFDNF AQTCANCHTQ DKAALQKVVA ERKQSINDLK IKVEDQLVHA HFEAKAALDA GATEAEMKPI QDDIRHAQWR WDLAIASHGI HMHAPEEGLR MLGTAMDKAA DARTKLARLL ATKGITHEIQ IPDISTKEKA QQAIGLNMEQ IKA EKQDFIK TVIPQWEEQA RKNGLLSQ
Source	E.coli
Target Names	nrfA
Protein Names	Recommended name: Cytochrome c-552 EC= 1.7.2.2 Alternative name(s): Ammonia-forming cytochrome c nitrite reductase Short name= Cytochrome c nitrite reductase
Expression Region	27-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.