



# Recombinant Escherichia coli O6:K15:H31 Cytochrome c-552 (nrfA)

<b>Product Code</b>	CSB-YP606953EGY
<b>Abbreviation</b>	nrfA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0T9X8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	EQTA APATPVTVEA KNETFAPQHP DQYLSWKATS EQSERVDALA EDPRLVILWA GYPFSRDYNK PRGHAFVTD VRETLRTGAP KNAEDGPLPM ACWSCSPDV ARLIQKDGED GYFHGKWARG GPEIVNNLGC ADCHNTASPE FAKGKPELTL SRPYAARAME AIGKPFKAG RFDQQSMVCG QCHVEYYFDG KNKAVKFPWD DGMKVENMEQ YYDKIAFSDW TNSLSKTPML KAQHPEYETW TAGIHGKNNV TCIDCHMPKV QNAEGKLYTD HKIGNPFDNF AQTCANCHTQ DKAALQKVVA ERKQSINDLK IKVEDQLVHA HFEAKAALDA GATEAEMKPI QDDIRHAQWR WDLAIASHGI HMHAPEEGLR MLGTAMDKAA DARTKLARLL ATKGITHEIQ IPDISTKEKA QQAIGLNMEQ IKAEEKQDFIK TVIPQWEEQA RKNGLLSQ
<b>Source</b>	Yeast
<b>Target Names</b>	nrfA
<b>Protein Names</b>	Recommended name: Cytochrome c-552 EC= 1.7.2.2 Alternative name(s): Ammonia-forming cytochrome c nitrite reductase Short name= Cytochrome c nitrite reductase
<b>Expression Region</b>	27-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.