



Recombinant Escherichia coli O6:K15:H31 Maltoporin (IamB)

Product Code	CSB-EP606954EGY
Abbreviation	IamB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TA25
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	VDFHG YARSGIGWTG SGGEQQCFQT TGAQSKYRLG NECETYAELK LGQEVWKEGD KSFYFDTNVA YSVAQKNDWE ATDPAFREAN VQGKNLIEWL PGSTIWAGKR FYQRHDVHMI DFYYWDISGP GAGLENIDVG FGKLSLAATR SSEAGGSSSF ASNNIYDYTN ETANDVFDVR LAQMEINPGG TLELGVDYGR ANLRDNYRLV DGASKDGWLF TAEHTQSVLK GFNKFVVQYA TDSMTSQGKG LSQGSGVAFD NEKFAYNINN NGHMLRILDH GAISMGDNWD MMYVGMYYQDI NWDNDNGTKW WTVGIRPMYK WTPIMSTVME IGYDNVESQR TGDKNNQYKI TLAQWQAGD SIWSRPAIRV FATYAKWDEK WGYDYTGSSS TNPYYGKAWS ADFNGGSFGR GDSDEWTFGA QMEIWW
Source	E.coli
Target Names	IamB
Protein Names	Recommended name: Maltoporin Alternative name(s): Maltose-inducible porin
Expression Region	26-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.