



# Recombinant Escherichia coli O6:K15:H31 Adenine deaminase (ade)

<b>Product Code</b>	CSB-MP606967EGY
<b>Abbreviation</b>	ade
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TB37
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNNSINHKFH HISRAEYQEL LAVSRGDAVA DYIIDNVSIL DLINGGEISG PIVIKGRYIA GVGAEYADAP ALQRIDARGA TAVPGFIDAH LHISSMMTP VTFETATLPR GLTTVICDPH EIVNVMGEAG FAWFARCAEQ ARQNQYLQVS SCVPALEGCD VNGASFTLEQ MLAWRDHPQV TGLAEMMDYP GVISGQNALL DKLDAFRHLT LDGHCPGLGG KELNAYIAAG IENCHESYQL EEGRRKLQLG MSLMIREGSA ARNLNALAPL INEFNSPQCM LCTDDRNPE IAHEGHIDAL IRRLIEQHNV PLHVAYRVAS WSTARHFGLN HLGLLAPGKQ ADIVLLSDAR KVTVQQVLVK GEPIDAQTLQ AEEARLAQS APPYGNTIDR QPVSASFAL QFTPGKRYRV IEVIHNELIT HSRSSVYSEN GFDRDDVCFI AVLERYGQRL APACGLLGGF GLNEGALAAT VSHDSHNIVV IGRSAEEMAL AVNQVIQDGG GLCVVRNGQV QSHLPLPIAG LMSTDTAQSL AEQIDALCAA ARECGPLPDE PFIQMAFLSL PVIPALKLTS QGLFDGEKFA FTTLEVTE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ade
<b>Protein Names</b>	Recommended name: Adenine deaminase Short name= Adenase Short name= Adenine aminase EC= 3.5.4.2
<b>Expression Region</b>	1-588
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.