



Recombinant Escherichia coli O6:K15:H31 Malate dehydrogenase (mdh)

Product Code	CSB-MP606981EGY
Abbreviation	mdh
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TCN0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKVAVLGAAG GIGQALALLL KTQLPSGSEL SLYDIAPVTP GVAVDLSHIP TAVKIKGFSG EDATPALEGA DVVLISAGVA RKP GMDRSDL FNVNAGIVKN LVQQVAKTCP KACIGIITNP VNTTVAIAAE VLKKAGVYDK NKLFGVTTLD IIRSNTFVAE LKGKQPGEVE VPVIGGHSGV TILPLLSQVP GVSFTEQEVA DLTKRIQNAG TEVVEAKAGG GSATLSMGQA AARFGLSLVR ALQGEQGVIE CAYVEGDGQY ARFFSQPLLL GKNGVEERKS IGTLSAFEQN ALEGMLDTLK KDIALGEEFV NK
Source	Mammalian cell
Target Names	mdh
Protein Names	Recommended name: Malate dehydrogenase EC= 1.1.1.37
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.