



Recombinant Escherichia coli O6:K15:H31 Ribosomal RNA large subunit methyltransferase G (rlmG)

Product Code	CSB-EP606984EGY
Abbreviation	rlmG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TD22
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSHLDNGFRS LTLQRFPATD DVNPLQAWEA ADEYLLQQLD DTEIRGPVLI LNDAFGALSC ALAEHKPYSI GDSYISELAT RENLRLNGID ESSVKFLDST ADYPQQPGVV LIKVPKTLAL LEQQLRALRK VVTPDTRIIA GAKARDIHTS TLELFEKVLG PTTTTLAWKK ARLINCTFNE PPLVDAPQTV SWKLEGTDWT IHNHANVFSR TGLDIGARFF MQHLPENLEG EIVDLGCGNG VIGLTLDDKN PQAKVVFVDE SPM AVASSRL NVETNMPEAL DRCEFMINNA LSGVEPFRFN AVLCNPPFHQ QHALTDNVAW EMFHHARRCL KINGELYIVA NRHLDYFHKL KKIFGNCTTI ATNNKFVVLK AVKLGRRR
Source	E.coli
Target Names	rlmG
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase G EC=2.1.1.174 Alternative name(s): 23S rRNA m2G1835 methyltransferase rRNA (guanine-N(2)-)-methyltransferase RlmG
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.