



Recombinant Escherichia coli O6:K15:H31 Nitric oxide reductase FIRd-NAD (+) reductase (norW)

Product Code	CSB-BP606993EGY
Abbreviation	norW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TEG9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSNGIVIIGS GFAARQLVKN IRKQDASIPL TLIAADSMDE YNKPDLSHVI SQGQRADDLT RQTAGEFAEQ FNLRLFPHTW VTDIDAEAHV VKSQNNQWQY DKLVLATGAS AFVPPVPGRE LILTLNSQQE YRACETQLRD ARRVLIVGGG LIGSELAMDF CRAGKAVTLI DNAASILASL MPPEVSSRLQ HRLTEMGVHL LLKSQLQGLE KTDSGILATL DRQRCIEVDA VIAATGLRPE TALARRAGLT INRGVCVDSY LQTSNADIYA LGDCAEINGQ VLPFLQPIQL SAMVLAKNLL GNNTPLKLPA MLVKIKTPEL PLHLAGETQR QDLRWHIYTE RQGMVARGVD DADQLRAFVV SEDRMKEAFG LLKTLSM
Source	Baculovirus
Target Names	norW
Protein Names	Recommended name: Nitric oxide reductase FIRd-NAD(+) reductase EC= 1.18.1.- Alternative name(s): Flavorubredoxin reductase Short name= FIRd-reductase Short name= FlavoRb reductase
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.